



Start of Transcription

Polyhedrin promoter

3901 TAGATCATGG AGATAATTAA AATGATAACC ATCTCGCAA TAAATAAGTA
wild-type ATG mutated to ATT

3951 TTTTACTGTT TCGTAACAG TTTTGTAATA AAAAAACCTA TAAATATTCC

Bam HI Rsr II BssH II

4001 GGATTATTCA TACCGTCCA CCATCGGGCG CGGATCCCGG TCCGAAGCGC

EcoRI Stu I Sal I Sst I Spe I Not I Nsp V

4051 GCGGAATTCA AAGGCCTACG TCGACGAGCT CACTAGTCGC GGCCGCTTTC

Xba I Pst I Xho I Sph I Kpn I Hind III

4101 GAATCTAGAG CCTGCAGTCT CGAGGCATGC GGTACCAAGC TTGTCGAGAA

SV40 polyadenylation signal

4151 GTACTAGAGG ATCATAATCA GCCATACCAC ATTGTAGAG GTTTTACTTG